



## OVENS AND INCUBATORS PREMIER RANGE

PREMIER  
SERIE

### MODELS:


- NATURAL AIR CONVECTION, DRYING AND STERILIZATION.
- FAN ASSISTED CIRCULATION, UNIVERSAL APPLICATIONS.
- NATURAL AIR CONVECTION, BACTERIOLOGY AND INCUBATION.

CONTROL: ANALOGUE OR DIGITAL MICROPROCESSOR CONTROL OF TEMPERATURE AND TIME, MODEL DEPENDENT.  
COMPLIES WITH THE STANDARDS: DIN 50011 - DIN 58945. REQUIRED FOR HEATING, STABILITY AND HOMOGENEITY.

### SAFETY:

STANDARD EN.61010. INCORPORATED FIXED OVER TEMPERATURE DEVICE .  
STANDARD DIN 12880. (CLASS 2 AND 3.1)SAFETY THERMOSTAT CONTROLLER FITTED.

### Leading edge technology



*Detailed longitudinal cross section.*

### COMMON FEATURES

#### Construction.

1. External case treated with a corrosive resistant epoxy coating.
2. Internal part: Easy to clean AISI 304 stainless steel double chamber, self adjusting door seal and adjustable shelves and guides.
3. Control panel: independent insulated control panel to facilitate all types of instruments, controls and regulators.
4. Adjustable air inlet.
5. Flexible floating door seal, self adjusting that maintains the best possible seal.

#### Technical Properties.

6. Excellent thermal qualities of the insulation has the optimum performance according to heater capacity and power consumption, with minimal external temperature loss.

7. Independent heating chamber for the heating elements to obtain an even heat distribution and rapid temperature equilibrium and stabilization.

Fan assisted convection models have a turbo fan. All incubators for bacteriology and cell culture have a second inner door of tempered glass.

#### Technology from J. P. Selecta:

8. Adjustable guide and shelf positions.

9. Double seal around the chamber to provide a gentle but effective seal.

10. Floating spring door that adjusts the pressure and absorbs the thermal expansion.

11. Adjustable door pressure system closure. Internal tempered glass door.


### NOTE:

For all models, the values for stability and homogeneity shown are based on temperature conditions with the ventilation closed. The optimum homogenization of temperature within the chamber is based on a reasonable load that does not surpass more than 70 % of the volume of the chamber. The graphic results shown for temperature for each model are based on the above criteria.

## CONTROL PANELS


### Models with Analogue control.

1. Main switch.
2. "On" indicator lamp.
3. Temperature control thermostat.
4. Heating "ON" indicator lamp.
5. Analogue thermometer temperature indicator.
6. Vacant positions for additional accessories.
7. Controllable safety thermostat that disconnects power to the heater in case of a fault in the main thermostat, manual reset (Directive DIN12880.2 class 2 and 3.1) and function signal lamp.



### Models with 4.3 inches TFT touch screen.

1. Main switch.
2. TFT touch screen:
  - Visual audible alarm .
  - Clock calendar.
  - Single or cyclic On / Off programming.
  - Up to 10 work programs.
  - Up to 6 segments per program.
  - Stability time in each segment (from 1 min to 99h) .
  - Alarms and events storage.
  - Probe error detection.
  - Self Diagnostics.
  - Ramps between segments.
  - Door open alarm.
  - Network failure detection and saving.
  - Over temperature and low temperature alarms and memorization (date, start time, end time and temperature).
  - Safety thermostat (TS) by software.
  - Mechanic safety thermostat (TS) .
  - PC software.
  - User manual on screen.
  - Temperature control auto-tuning.
  - Configurable parameters: Date / time, temperature correction , data collection interval, language (English, Spanish and French) , °C / °F selection , over temperature and low temperature limit.
3. RS-232 output.
4. USB output.
5. Security thermostat.
6. Ethernet output para for LAN connection.



## MODEL SUMMARY TABLE

| Models                 | CONTERM                 | DIGITHEAT               | DIGITRONIC         | INCUBAT                   | INCUDIGIT                 |
|------------------------|-------------------------|-------------------------|--------------------|---------------------------|---------------------------|
| <b>TYPE</b>            | Drying Oven             | Drying Oven             | Universal          | Bacteriological Incubator | Bacteriological Incubator |
| <b>CONTROL</b>         | Temperature             | Temperature + time      | Temperature + time | Temperature               | Temperature + time        |
| <b>DISPLAY</b>         | Analogue                | Digital                 | Digital            | Analogue                  | Digital                   |
| <b>AIR</b>             | Convection              | Convection              | Fan assisted       | Convection                | Convection                |
| <b>CIRCULATION</b>     | natural                 | natural                 |                    | natural                   | natural                   |
| <b>CAPACITY LITRES</b> | 19 - 36 - 52 - 80 - 150 | 19 - 36 - 52 - 80 - 150 | 33 - 47 - 76 - 145 | 19 - 36 - 52 - 80 - 150   | 19 - 36 - 52 - 80 - 150   |

## ACCESSORIES



Part No.  
**2000002** Timer switch 0-120 minutes.  
 Suitable for **CONTERM**.  
**2000003** Timer switch 0-12 hours.  
 Suitable for **CONTERM** and **INCUBAT**.

**2000009** 24 hour programmer with continuous on/off cycling up to every 15 minutes.  
 Suitable for **CONTERM** and **INCUBAT**.



Part No.  
**2000016** Digital printer for time and temperature with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.  
 Suitable for **DIGITHEAT**, **DIGITRONIC** and **INCUDIGIT**.



**Optional communication modules**  
 Part No. **2101623** Module for Wifi network.  
 Part No. **2101624** Module for Bluetooth.  
 Part No. **2101625** Module RF.  
 Part No. **2101626** RS-232 to RS-485 converter.  
 Suitable for **DIGITHEAT**, **DIGITRONIC** and **INCUDIGIT**.



## Incubators for bacteriology and cell culture “Incubig-TFT”

**NATURAL CONVECTION.**  
**MICROPROCESSOR CONTROL AND DIGITAL DISPLAY OF TEMPERATURE AND TIME.**  
**ADJUSTABLE TEMPERATURE FROM AMBIENT +5 °C TO 80 °C.**  
**STABILITY: ±0.1 °C, UP TO 37 °C. HOMOGENEITY: ±0.5 °C, UP TO 37 °C.**  
**SET ERROR: ±2% OF THE WORKING TEMPERATURE, RESOLUTION 0.1 °C**  
**INTERNAL TEMPERED GLASS DOOR.**

PREMIER  
SERIE



NEW  
DESIGN



### SAFETY:

STANDARD EN.61010 OVER TEMPERATURE CUT OUT FITTED.  
 STANDARD DIN 12880. ADJUSTABLE SAFETY THERMOSTAT FITTED.

### Capacities up to 720 litres

#### FEATURE


Microprocessor control and 4.3 inches TFT touch screen display.  
 Large surface area heating elements.  
 Inner chamber made of AISI 304 stainless steel.  
 Double door, interior door of tempered glass that allows the user to see the contents of the chamber without opening the door.  
 Adjustable air vent.  
 Epoxy covered external case.

#### STANDARD EQUIPMENT

For Part No. 2000238, 2 shelves and 4 shelf guides.  
 For Part No. 2000239 and 2000240, 2 shelves.



Model Part No. 2000238.



Performance graph of temperature and time.

- A. Set at 80 °C: 1 h 45'.
- B. Set at 56 °C: 1 h 10'.
- C. Set at 37 °C: 54'.

Note: To obtain the optimum homogeneity at the set temperature, the load should not surpass more than 70 % of the volume of the chamber.



Models Part No. 2000239 and 2000240.

## CONTROL PANEL

### 4.3 inches TFT touch screen models:

1. Main switch.
  2. TFT touch screen:
- Visual audible alarm.
  - Clock calendar.
  - Single or cyclic On / Off programming.
  - Up to 10 work programs.
  - Up to 6 segments per program.
  - Stability time in each segment (from 1 min to 99h).
  - Alarms and events storage.
  - Probe error detection.
  - Self Diagnostics.
  - Ramps between segments.
  - Door open alarm.
  - Network failure detection and saving.
  - Over temperature and low temperature alarms and

memorization (date, start time, end time and temperature).

- Safety thermostat (TS) by software.
  - Mechanic safety thermostat (TS).
  - USB and RS -232 output.
  - PC software.
  - User manual on screen.
  - Temperature control auto-tuning.
  - Configurable parameters: Date / time, temperature correction, data collection interval, language (English, Spanish and French), °C / °F selection, over temperature and low temperature limit.
3. RS-232 output.
  4. USB output.
  5. Security thermostat.
  6. Ethernet output para for LAN connection.



## MODELS

| Part No. | Type   | Capacity litres | Height / Width / Depth (interior) cm |     |    | Height / Width / Depth (exterior) cm |     |    | N° of shelf guides | Power W | Weight Kg |
|----------|--------|-----------------|--------------------------------------|-----|----|--------------------------------------|-----|----|--------------------|---------|-----------|
| 2000238  | 1 door | 288             | 80                                   | 60  | 60 | 97                                   | 91  | 76 | 8                  | 570     | 87        |
| 2000239  | 2 door | 400             | 100                                  | 80  | 50 | 130                                  | 114 | 75 | 10                 | 1100    | 160       |
| 2000240  | 2 door | 720             | 120                                  | 100 | 60 | 152                                  | 134 | 85 | 12                 | 1600    | 225       |

## ACCESSORY



**4120131 USB adapter model.**  
Pen-Drive included (Memory board) for data storage.

Accessories must be factory installed.



**2000016 Digital printer for time and temperature** with numerical printout on continuous paper roll, with print intervals from 1 minute to 99 hours.

## SPARES

Shelves and guides.

| Oven Part No.           | 2000238        | 2000239        | 2000240        |
|-------------------------|----------------|----------------|----------------|
| <b>Shelves</b>          | <b>2002372</b> | <b>2000063</b> | <b>2000064</b> |
| <b>Guides (2) (Set)</b> | <b>2002371</b> | -              | -              |

Each self requires two guides i.e. one set.



## Incubator for Petri capsules

NATURAL CONVECTION.

MICROPROCESSOR REGULATION AND TEMPERATURE DIGITAL CONTROL.

FOR ADJUSTABLE TEMPERATURES FROM AMBIENT +5°C TO 60°C.

STABILITY: ±0,1°C TO 37°C. HOMOGENEITY: ±0,1°C TO 37°C. SETPOINT ERROR: ±0,1°C. RESOLUTION: 0,1°C.

**Small size. Culture visual control. Transportable.**

## APPLICATIONS

Specially designed for bacteria and fungi cultures in Petri capsules at the same temperature of human body.

## FEATURES

Culture surface 320 x 220 mm (Inner height: 20mm)

Culture visual monitoring.

Transparent cover.

Easy access to samples.

**Approximate capacity:** (single level) (mm)

15 Petri capsules of Ø55.

10 Petri capsules of Ø80.

7 Petri capsules of Ø90.

6 Petri capsules of Ø100.

3 Petri capsules of 120x120.

2 Petri capsules of Ø140.

## MODEL

| Part No. | Margin T° °C | Height / Width / Depth (Exterior) cm |    |    | Height / Width / Depth (interior) cm |    |    | Power W. | Weight Kg |
|----------|--------------|--------------------------------------|----|----|--------------------------------------|----|----|----------|-----------|
| 4002629  | Amb +5 to 60 | 7                                    | 45 | 25 | 2                                    | 32 | 22 | 300      | 3.5       |

